
# Percussion Massage Belt™

User's Manual / Manual de Usuario / Guide d'Utilisation

# Percussion Massage Belt™ by Wagan Tech®

---

Thank you for purchasing the Percussion Massage Belt by Wagan Tech. With normal care and proper treatment, it will provide years of reliable service. Please read all operating instructions carefully before use.

## FEATURES

- 2 Rhythmic Massagers
- 10 Intensity levels
- 20 Different Combinations
- Versatile: Use for Neck/Shoulders, Back, Lumbar, Legs and Posterior

## SPECIFICATIONS

Rated Voltage	110V
Frequency	60Hz
Power	60W

## WARNINGS

**To reduce the risk of burning, fire, electric shock, or physical injury, please follow the precautions below.**

- Do not unplug the unit when it is still operating. Power off the unit before unplugging.
- Always unplug the unit when it is not in use and when there is a power outage.
- Do not carry the unit by the power cord.
- Never operate the unit with a damaged power cord or plug. Do not attempt to repair the power cord if it is damaged. Please contact us for assistance.
- Do not spill liquid on any part of the unit. Do not use the unit in the bathroom or other damp and humid places. This may cause electric shock or cause the unit to malfunction.

- Keep the unit and power cord away from heat sources, including electric blankets.
- Do not use the unit with a transformer as it may lead to malfunction or electric shock.
- To avoid electric shock or injury, do not attempt to open or disassemble any part of the unit.

### PERSONAL SAFETY

- Do not stand or sit on the unit. Only use the unit for the purposes described in this manual.
- Do not use the unit against bare skin: while thin clothing may increase effectiveness, exposing the skin directly to the massage may irritate the skin.
- Persons with any of the conditions listed below or undergoing medical treatment should consult their physician before using the unit:
  - » Pregnant women, people who are ill, or those in poor physical shape.
  - » People who are prohibited from receiving a massage by a physician due to thrombosis, aneurysm, acute varicose veins, or other circulatory disorders.
  - » People with irregular curvatures of the spine.
- If you start feeling sick or if the massage seems painful during operation, please discontinue use immediately.
- Do not fall asleep while using the unit.
- Do not use the unit after consuming alcohol.
- Do not allow children to play with the unit.

## OPERATION

The Percussion Massage Belt has a safety feature causing the massage motors to shut off if overheated. If one or both massage motors does not activate, allow unit to cool completely before attempting to use again.


1. Plug the power cord into a standard wall outlet. After a beep, the appliance will go into standby mode and the POWER LED will blink red.
2. Place the Percussion Massage Belt on your shoulders, back, or legs as shown in the Massage Areas section.
3. *Automatic Mode:*  
Press the **(M)** (mode) button and the unit will switch to automatic mode. The "AUTO" indicator will light green. The massager will automatically cycle through all 20 different massage rhythm combinations over a 15 minute program.

### *Manual Mode:*


- Press the **(M)** button a second time and the unit will switch to manual mode. The "MANUAL" indicator will light red. The massager will run a single massage combination continuously for 15 minutes. You can select a different massage combination by pressing the **(M)** button repeatedly until you have reached the desired pattern.
4. There are 10 levels massaging intensity. Press the **(+)** and **(-)** buttons anytime during operation for adjustment.
  5. To stop the unit before its cycle is over, press and hold the **(M)** button for 3 seconds, and the unit will turn off. For both Automatic and Manual modes, the massager's timer will automatically shut off the unit after its 15-minute cycle ends. Allow the unit to cool for at least 15 minutes before reusing.

## MASSAGE AREAS


The unit massages the shoulder, back, lumbar, legs, or posterior with rhythmic percussion. See illustrations below for usage examples.


Neck/Shoulders


Upper Back

Lumbar (Lower Back)

Legs

Posterior

### MAINTENANCE

- Always unplug the unit before cleaning it. Never touch the power plug with wet hands.
- Wipe the unit with a soft, dry cloth. Do not use cloths containing any kind of chemical or alcohol and solvent liquids.
- Never submerge any part of the unit in liquid.
- If the unit is particularly dirty, use the following instructions:
  1. Dampen a cloth in water or a 3%–5% mild detergent solution.
  2. Wipe the dirty areas with the wet cloth.
  3. Clean the cloth in water and wring it out, then wipe away any remaining detergent off the surface of the unit.
  4. Wipe the areas with a soft, dry cloth to dry.
  5. Wait for the unit to become fully dry before use.
- Do not press the unit against plastic for a long time.

# WAGAN Corp. Limited Warranty

**The WAGAN Corporation warranty is limited to products sold only in the United States.**

**Warranty Duration:**

Product is warranted to the original purchaser for a period of one (1) year from the original purchase date, to be free of defects in material and workmanship. WAGAN Corporation disclaims any liability for consequential damages. In no event will WAGAN Corporation be responsible for any amount of damages beyond the amount paid for the product at retail.

**Warranty Performance:**

During the warranty period, a product with a defect will be replaced with a comparable model when the product is returned to WAGAN Corporation with an original store receipt. WAGAN Corporation will, at its discretion, replace or repair the defective part. The replacement product will be warranted for the balance of the original warranty period. This warranty does not extend to any units which have been used in violation of written instructions furnished.

**Warranty Disclaimers:**

This warranty is in lieu of all warranties expressed or implied and no representative or person is authorized to assume any other liability in connection with the sale of our products. There shall be no claims for defects or failure of performance or product failure under any theory of tort, contract or commercial law including, but not limited to negligence, gross negligence, strict liability, breach of warranty, and breach of contract.

**Returns:**

Please contact customer service by phone, email, or visiting our website under "Customer Care" page for instructions on how to process a warranty claim. WAGAN Corporation is not responsible for any item(s) returned without an official Return Authorization number (RA#).

Register your product online at <http://tinyurl.com/wagan-registration> to be added to our email list. You will receive previews on our upcoming products, promotions, and events.

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©2015 Wagan Corporation. All Rights Reserved
Wagan Tech and *wagan.com* are trademarks used by Wagan Corporation
# Garantie Limitée de WAGAN Corporation
## **La garantie de WAGAN Corporation est limitée aux produits vendus uniquement aux États-Unis.**
Hayward, CA 94544, U.S.A.
© 2015 Wagan Corporation. Tous droits réservés.